

## Product Specification

---

|               |                            |
|---------------|----------------------------|
| Product Name: | Glutathione oxidized       |
| Grade:        | High Purity                |
| Product Code: | GB0228                     |
| Formula:      | $C_{20}H_{32}N_6O_{12}S_2$ |
| MW:           | 612.63                     |
| CAS#:         | 27025-41-8                 |

---

---

| Test Items                        | Specifications           |
|-----------------------------------|--------------------------|
| Appearance                        | White Crystalline Powder |
| Specific Rotation $[\alpha]^{20}$ | -94° - -102°             |
| Recommended Retest Period:        | 2 Years                  |

---

Storage Condition: 2-8°C. Keep dry. Warm to room temperature before opening.

*Rev 0*

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)