
Product Specification

Product Name: Beef Extract
Product Code: BB0114

Product Description:

Beef Extract: This product is prepared from very low fat bovine skeletal muscle and marrow, which is spinal column free with no tendons. Beef Extract Powder is highly nutritious and is used in preparing microbiological culture media. It is recommended for use in culture media for bacteriological examination of water, milk and other materials.

Test Items	Specifications
Amino Nitrogen (AN)	≥ 3.9 %
Total Nitrogen (TN)	≥ 10.0 %
Loss on drying	≤ 6 %
Ash	≤ 16 %
pH (2% solution,@25°C)	6.5 – 7.5
Standard plate count	< 5000 col/g
Yeasts and molds	<100 col/g
Coliforms	Negative
<i>Salmonella</i>	Negative

Storage Condition: 18 – 25°C, Protect from moisture