
Product Specification

Product Name: LB Broth Powder, Lennox
Product Code: SD7005

Product Description:

LB Broth Powder, Lennox: Lennox LB is a highly-referenced microbial growth medium used for the cultivation of E. coli. This nutrient-rich microbial broth contains peptides, amino acids and carbohydrates in a low-salt formulation.

Components:

Tryptone – 10g
Yeast Extract – 5g
Sodium Chloride – 5g

Test Items

Specifications

Physical and Chemical Tests

DCM Appearance	Homogeneous, free flowing, light beige
pH (2.5% solution, @25°C)	6.00 – 7.50

Storage Condition: 18 – 25°C