

---

## Product Specification

---

Product Name: Peptone-A, (from Meat)  
Product Code: G213

---

Product Description:

A specialized peptone prepared by the enzymatic digestion of selected fresh meat to be used in culture media.

---

Test Items	Specifications
Amino Nitrogen (AN)	≥2.6%
Total Nitrogen (TN)	≥12.0%
Loss on drying	≤6.0%
Ash	≤15.0%
pH (2% solution,@25°C)	6.5-7.5
Standard Plate count	<5000 CFU/g
Molds and Yeast	<100 CFU/g
Coliforms	Negative
<i>Salmonella</i>	Negative
<u>Growth Supporting Properties</u>	
Peptone Agar(USP)	Satisfactory
Broth Formulation	Satisfactory

---

Storage Condition: 18 – 25°C, Protect from Moisture