
Product Specification

Product Name: Polypeptone
Product Code: G258

Product Description:

PolyPeptone: This peptone is a mixture of Casein Peptone and Meat Peptone. It can be used in general culture media as a superior source of nutrients. The growth of some microorganisms may be better than when the individual peptones are used.

Test Items

Specifications

Appearance (dry)	Fine, homogenous, light to medium tan yellow
Total Nitrogen (TN)	≥10.0%
Loss on drying	≤6.0%
Ash	≤15.0%
pH (2% solution, @25°C)	6.5-7.5
Standard Plate count	<5000 CFU/g
Molds and Yeast	<100 CFU/g
Coliforms	Negative
<i>Salmonella</i>	Negative

Storage Condition: 18 – 25°C, Protect from Moisture