
Product Specification

Product Name: SSC Buffer – 20X Liquid Concentrate
Grade: Biotechnology
Product Code: D0623

Product Description:

SSC Buffer is a standard reagent used in Southern and Northern hybridization procedures to control the stringency of the wash buffer upon usage in the washing steps after hybridization.

Note: This product (20X) consists of 3.0M Sodium Chloride and 0.3M Sodium Citrate (adjusted to pH 7.0 with HCl).

Test Items	Specifications
------------	----------------

Appearance	Clear Solution
pH @25°C	6.85 - 7.15

Recommended Retest Period: 1 Year

Storage Condition: 18 - 25°C, Protect from Moisture

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com