
Product Specification

Product Name: TBS Buffer (Tris Buffered Saline), 10X Solution
Grade: Ultra Pure, 0.22 μ M filtered
Product Code: UA00279

Product Description:

This product is ideal for use as a wash buffer and antibody diluent for ELISA, Western blotting and other immunoassays.

Note: The salt concentration in this product in this buffer is at 10X. A single (1X) strength is 25mM Tris-HCl, 150mM NaCl.

Test Items	Specifications
------------	----------------

Appearance	Clear Solution
pH (1X, @25°C)	7.2 – 7.6
Sterility	Passes
DNase	None
RNase	None
Protease	None

Recommended Retest Period: 1 Year

Storage Condition: 18 - 25°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com