
Product Specification

Product Name: TBST (TBS-Tween 20) 10X
Grade: Ultra Pure, 0.22 μ M filtered
Product Code: UA0046

Product Description:

This product is commonly used as wash buffer in Western Blot, ELISA, and immunoassays. It can also be used as blocking buffer for plate based assays.

Note: The salt concentration of this buffer is at 10X. A single (1X) strength contains 0.025M Tris, 0.15M NaCl, 0.1% Tween-20.

Test Items	Specifications
------------	----------------

Appearance	Clear Solution
pH (1X, @25°C)	7.4 – 7.8
Sterility	Passes
DNase	None
RNase	None
Protease	None

Recommended Retest Period: 1 Year

Storage Condition: 18 - 25°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com