
Product Specification

Product Name: TG Buffer (Tris-Glycine) 10X Solution, pH 8.3
Grade: Biotechnology
Product Code: A0028

Product Description:

This product is used to make a transfer buffer in Western blotting and gel electrophoresis.

Note: The salt concentration of the buffer is at 10X. A single strength (1X) solution contains 0.025M Tris base and 0.192M Glycine.

Test Items	Specifications
Appearance pH (1X, @25°C)	Clear Solution 8.10 – 8.50

Recommended Retest Period: 1 Year

Storage Condition: 18 - 25°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com