

---

## Product Specification

---

Product Name:	Agarose A
Grade:	Biotechnology
Product Code:	D0012
CAS#:	9012-36-6

---

---

Test Items	Specifications
------------	----------------

---

Appearance	Fine, white homogeneous powder
Moisture content	≤10%
Gel Strength 1.5%	≥1,120 g/cm <sup>2</sup>
Gelling Point	34.5-37.5°C
Electroendosmosis-Mr	0.09 - 0.13
Sulfate	≤0.15%
DNase, RNase and Protease	None detected

Recommended Retest Period:	2 Years
----------------------------	---------

---

Storage Condition: 18 – 25°C, Protect from moisture

*Rev 0*

---

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)