

Product Specification

Product Name:	Albumin, bovine, fraction V, heat shock isolation pH 7.0
Grade:	Biotechnology
Product Code:	AD0023
Formula:	N/A
MW:	N/A
CAS#:	9048-46-8

Test Items	Specifications
Appearance	Off-white to pale green/yellow lyophilized powder
Purity	≥98.0%
Protein content	≥97.0%
Residue on Ignition	≤2.0%
pH (10%, H ₂ O)	6.5 - 7.5
Loss on drying	≤5.0%
Solubility(@4% in H ₂ O)	Freely soluble

Note: DNase, RNase & Protease are not tested

Storage Condition: 2 – 8°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com