

Product Specification

Product Name:	EDTA, Disodium Salt, Dehydrate
Grade:	Ultra Pure
Product Code:	EB0185
Formula:	$C_{10}H_{14}N_2O_8Na_2 \cdot 2H_2O$
MW:	372.24
CAS#:	6381-92-6

Test Items

Specifications

Appearance	White crystalline powder
Assay	≥99.0%
Reaction of water solution	Passes
Clarification test	Passes
pH (1%)	4.0 - 5.0
pH (5%)	4.0 - 6.0
Chloride	≤0.01%

Recommended Retest Period: 2 Years

Storage Condition: 18 – 25°C, Protect from moisture

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com