

Product Specification

Product Name:	EGTA
Grade:	High Purity
Product Code:	ED0077
Formula:	$C_{14}H_{24}N_2O_{10}$
MW:	380.35
CAS#:	67-42-5

Test Items	Specifications
Appearance	White crystalline powder
Assay	99.0 - 101.0%
Heavy metals	≤ 5ppm
Melting point	240 - 244°C
Loss on drying	≤ 2.0%
Chloride	≤ 50ppm
Recommended Retest Period:	2 Years

Storage Condition: 18 – 25°C, Protect from moisture

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com