

---

## Product Specification

---

Product Name:	NP-40 (octylphenoxypolyethoxylethanol)
Grade:	Reagent
Product Code:	NDB0385
Formula:	N/A
MW:	603
CAS#:	9016-45-9

---

---

### Test Items

### Specifications

---

Appearance	Clear to Hazy Colorless to Slight Yellow Liquid
Solubility (5% w/v in H <sub>2</sub> O)	Passes
Density at 20°C (g/mL)	1.02 - 1.10
Cloud point (11%, aq.)	52 - 56°C
Recommended Retest Period:	2 Years

---

Storage Condition: 18 – 25°C

*Rev 0*

---

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)