
Product Specification

Product Name:	Sodium Chloride
Grade:	Biotechnology
Product Code:	DB0483
Formula:	NaCl
MW:	58.44
CAS#:	7647-14-5

Test Items	Specifications
Appearance	White Crystalline powder
Assay	≥99.5%
Iodide	≤0.002%
Phosphate	≤0.0005%
Calcium	≤0.005%
Iron	≤0.0002%
Heavy Metals	≤0.0005%
Recommended Retest Period:	2 Years

Storage Condition: 18 – 25°C, Protect from moisture

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com