

---

## Product Specification

---

Product Name:	Streptavidin
Grade:	Biotechnology
Product Code:	SE497
Formula:	N/A
MW:	N/A
CAS#:	9013-20-1

---

---

Test Items	Specifications
------------	----------------

---

Appearance	White Freeze-dried Powder
Protein Content	650 ± 150 mg/g powder
Bioactivity	≥17Units/mg
SDS-PAGE Electrophoresis Purity	≥95%
HPLC Purity	≥95%
pH	6.5 ± 0.5
Molecular Weight	53.0 ± 5.3 KDa
Bacterial Endotoxins	≤1.0EU/ug rSTA

---

Storage Condition: -20°C, Protect from moisture. Warm to room temperature before opening.

*Rev 0*

---

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)