

## Product Specification

---

Product Name:	TCEP hydrochloride
Grade:	High Purity
Product Code:	TB0974
Formula:	$C_9H_{15}O_6P.HCl$
MW:	286.6
CAS#:	51805-45-9

---

---

Test Items	Specifications
Appearance	White crystals
Infrared Spectrum	Conforms reference spectrum
Solubility 30%	Clear colorless aqueous solution
TLC	≥99%
Functional reduction(DTNB)	≥95%
Recommended Retest Period:	2 Years

---

Storage Condition: 2 – 8°C

*Rev 0*

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)