
Product Specification

Product Name:	Triton X-100
Grade:	Biotechnology
Product Code:	TB0198
Formula:	$C_{14}H_{22}O(C_2H_4O)_n(n=9 - 10)$
MW:	647 (Avg.)
CAS#:	9002-93-1

Test Items	Specifications
Appearance	Hazy viscous liquid
Cloud Point (1% aq.)	63-69
Color (APHA)	≤60
Recommended Retest Period:	2 Years

Storage Condition: 18 – 25°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com