

---

## Product Specification

---

Product:	dNTP Mixture, 25mM
Grade:	Biotechnology
Formula:	N/A
MW:	N/A
CAS#:	N/A
Product Code:	DD0057

---

---

Test Items	Specifications
------------	----------------

---

Appearance	Clear and colorless solution
Dnase, Rnase	None detected
pH	6.9 - 7.1

---

Description: A solution that is 25mM dATP, 25mM dGTP, 25mM dCTP, 25mM dTTP.

Storage: Long term (infrequent use, 1-2 times per month): -70°C.

Daily/Weekly use: -20°C. Always avoid freeze-thaw cycles or exposure to frequent temperature changes. These fluctuations can greatly alter product stability.

---

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: [sg@biobasic-asia.com](mailto:sg@biobasic-asia.com) Web: [www.biobasic-asia.com](http://www.biobasic-asia.com)