
Product Specification

Product Name: Gentamicin Sulfate Solution (15mg/ml, Sterile)
Product Code: BSAT006010

Product Description:

Gentamicin Sulfate Solution is an aminoglycoside antibiotic that can cause misreading of the mRNA sequence and inhibition of translocation, inhibiting protein synthesis. It is a broad-spectrum antibiotic that is effective against Gram-positive and Gram-negative bacteria. In tissue culture, it can be used to inhibit the growth of various strains of mycoplasma

Product Content(s)

Gentamicin Sulfate Solution (15mg/ml, Sterile)	10ml
---	------

Note: Store at -20°C. This reagent has a shelf-life of 2 years.

Things to take note:

1. Product is dissolved in distilled water
2. Product should be handled in aseptic conditions to avoid contamination
3. Product should not to be placed at room temperature for a long time or stored at 4°C for extended period
4. Thaw at 2°C to 8°C, shake well before use, and avoid repeated freezing and thawing. It is recommended to aliquot the reagent into smaller volume and keep it in the freezer

Note:

This product is designed for research use only. Not to be used on humans or animals' diagnosis and other in vivo experiments. Please wear safety gear when using.

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com