

---

### Product Specification

---

Product Name: Neomycin Sulfate Solution (10mg/ml, Sterile)  
Product Code: BSAT004010

---

#### Product Description:

Neomycin Sulfate Solution is an aminoglycoside antibiotic that can cause misreading of the mRNA sequence and inhibition of translocation, inhibiting protein synthesis. It is a broad-spectrum antibiotic that is effective against Gram-positive and Gram-negative bacteria. It is used as an antimicrobial agent to prevent contamination, as well as a selection agent for prokaryotic cells that have been transformed using *neo* as a selectable marker gene.

---

#### Product Content(s)

---

Neomycin Sulfate Solution (10mg/ml, Sterile)	10ml
---	------

---

Note: Store at -20°C. This reagent has a shelf-life of 2 years.

#### Things to take note:

1. Product is dissolved in distilled water
2. Product should be handled in aseptic conditions to avoid contamination
3. Product should not to be placed at room temperature for a long time or stored at 4°C for extended period
4. Thaw at 2°C to 8°C, shake well before use, and avoid repeated freezing and thawing. It is recommended to aliquot the reagent into smaller volume and keep it in the freezer

#### Note:

This product is designed for research use only. Not to be used on humans or animals' diagnosis and other in vivo experiments. Please wear safety gear when using.