
Product Specification

Product Name:	Nitrocellulose Membrane, Pore Size 0.45 μ m
Product Code:	BEW00829
Dimensions	300mm x 3000mm
Size(s):	1 roll

Test Items	Specifications
BSA Binding Capacity	125 μ g/cm ²
Dimensions	300mm x 3000mm
Pore Size	0.45 μ m
Thickness	150 μ m \pm 10 μ m

Storage Condition: 18 – 25°C

Rev 0

2 International Business Park Road Strategy #01-04, Singapore(609930)

Tel: (+65) 6954 2519, (+65) 6491 5938

Email: sg@biobasic-asia.com Web: www.biobasic-asia.com