

---

## Product Specification

---

**Product Name:** 0.1% Brilliant Green Aqueous Solution

**Product Code:** SG7253

---

### PRODUCT DESCRIPTION:

Supplement for Tetrathionate Broth Base (TTB).

### USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Aseptically add the contents of one vial each of the 0.1% Brilliant Green Aqueous Solution and Iodine Solution to every 100 mL of TTB medium to prepare TTB enrichment broth.

---

---

Item	Specification
<b>Packing Size:</b>	1 mL x 10 Vials/box
<b>Formulation:</b>	1 mg Brilliant Green
<b>Storage Condition:</b>	2-8°C
<b>Shelf Life:</b>	1 Year

---

**Waste Disposal:** After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.