
Product Specification

Product Name: 50% Egg Yolk Suspension

Product Code: SG7353

PRODUCT DESCRIPTION:

Acts as a supplement for MYP Agar Base (Mannitol egg yolk polymyxin).

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial each of the 50% Egg Yolk Suspension and 1 mg Polymyxin B (SG7352) to every 95 mL of MYP Base medium to prepare MYP medium.

Item	Specification
Packing Size:	5 mL x 10 Vials/box
Formulation:	Physiological saline, egg yolk
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be treated by autoclaving at 121°C for 30 minutes.