
Product Specification

Product Name: Brucella Broth

Product Code: SG7380

PRODUCT DESCRIPTION:

Brucella Broth is used for the cultivation and isolation of Campylobacter.

Use:

Stir to suspend 28.1 g of dehydrated medium in 1L distilled H₂O. Mix well with heat with frequent agitation. Boil for one minute. Autoclave for 15 minutes at 121°C to sterilize. Allow to cool (40 – 45°C) before making additions, such as antibiotics (if desired), before dispensing. For laboratory use only.

Components:

Peptone – 10.0g/L

Yeast extract – 2.0g/L

Casein peptone – 10.0g/L

Sodium chloride – 5.0g/L

Glucose – 1.0g/L

Sodium bisulfite – 0.1g/L

Item	Description
------	-------------

Physical and Chemical Tests

Appearance: Yellowish Powder

pH (@25oC): 6.8 – 7.2

Storage: Store in a dry and dark place at room temperature. Tighten the cap immediately after use.