

---

## Product Specification

---

**Product Name:** Candida Chromogenic Agar

**Product Code:** SG7148

---

### PRODUCT DESCRIPTION:

Candida Chromogenic Agar is used for the qualitative direct detection, differentiation and presumptive identification of Candida species.

### Use:

Stir to suspend 47.7 g of dehydrated medium in 1L of distilled or deionized Water. Mix well with heat with frequent agitation. Boil to completely dissolve the powder, then cool to about 50°C. Do not overheat or autoclave. Pour into Petri dishes and set aside for use. For laboratory use only.

### Components:

Peptone – 10.2g/L

Agar – 15.0g/L

Chloramphenicol – 0.5g/L

Chromogenic mix – 22.0g/L

---

Item	Description
------	-------------

---

### Physical and Chemical Tests

**Appearance:** Yellowish Powder

**pH (@25°C):** 6.1±0.2

---

**Storage:** Store in a dry and dark place at room temperature. Tighten the cap immediately after use.