
Product Specification

Product Name: Cefixime

Product Code: SG7223

PRODUCT DESCRIPTION:

Supplement for the Modified MacConkey Sorbitol Agar.

USE:

The vial has an arrow mark on the aluminum cap. Before opening the vial, disinfect the surface of the vial with 75% alcohol. Under aseptic conditions, open the aluminum cap in the direction of the arrow, tear off the aluminum cap, and open the rubber stopper of the vial. Aseptically add the contents of one vial of Potassium Tellurite Solution (SG7224) and 5 vials of Cefixime to every 500 mL modified sorbitol MacConkey (CT-SMAC) agar medium. Alternatively, add one vial of Modified Sorbitol MacConkey Agar Additive (SG7222) to 200mL of the agar medium.

Item	Specification
Packing Size:	2 mL x 10 Vials/box
Formulation:	5 µg Cefixime
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.