
Product Specification

Product Name: Ceftazidime

Product Code: SG7285

PRODUCT DESCRIPTION:

Acts as a supplement for PALCAM Agar Base. Ceftazidime provides a broad spectrum of inhibition against other potential contaminants.

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial of the solution along with one vial each of 0.5 mg Acriflavine (SG7283) and 1 mg Polymyxin B (SG7284) to every 100 mL of PALCAM Agar Base Medium.

Item	Specification
Packing Size:	1 mL x 10 Vials/box
Formulation:	2.0 mg Ceftazidime
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.