
Product Specification

Product Name: Coliform Escherichia Coli (ECC) Chromogenic Agar

Product Code: SG7146

PRODUCT DESCRIPTION:

Coliforms & E. coli (ECC) Chromogenic Agar is used for the rapid cultivation of E. coli and Coliforms.

Use:

Stir to suspend 32.8 g of dehydrated medium in 1L of distilled or deionized Water. Mix well with heat with frequent agitation. Boil to completely dissolve the powder, then cool to about 50°C. Do not overheat or autoclave. Pour into Petri dishes and set aside for use. For laboratory use only.

Components:

Nutrient Components – 17.0g/L

Agar – 15.0g/L

Chromogenic Mixture – 0.8g/L

Item	Description
------	-------------

Physical and Chemical Tests

Appearance: Yellowish Powder

pH (@25°C): 7.0±0.2

Storage: Store in a dry and dark place at room temperature. Tighten the cap immediately after use.