
Product Specification

Product Name: D- Cycloserine

Product Code: SG7366

PRODUCT DESCRIPTION:

Acts as a supplement for Tryptose Sulfite Cycloserine (TSC) Agar Base.

USE:

Add one vial of D- Cycloserine to every 100 mL of TSC Agar Base.

Item	Specification
Packing Size:	5 mL x 10 Vials/box
Formulation:	0.04 g D-Cycloserine
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be treated by autoclaving at 121°C for 30 minutes.