
Product Specification

Product Name: Iodine Solution

Product Code: SG7252

PRODUCT DESCRIPTION:

Supplement for Tetrathionate Broth Base (TTB).

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Aseptically add the contents of one vial each of the 0.1% Brilliant Green Aqueous Solution and Iodine Solution to every 100 mL of TTB medium to prepare TTB enrichment broth.

Item	Specification
Packing Size:	2 mL x 10 Vials/box
Formulation:	0.5 g Potassium Iodide, 0.4 g Iodine, 2 mL Distilled Water
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.