
Product Specification

Product Name: Modified Sorbitol MacConkey Agar Additive

Product Code: SG7222

PRODUCT DESCRIPTION:

Supplement for the Modified MacConkey Sorbitol Agar.

USE:

The vial has an arrow mark on the aluminum cap. Before opening the vial, disinfect the surface of the vial with 75% alcohol. Under aseptic conditions, open the aluminum cap in the direction of the arrow, tear off the aluminum cap, and open the rubber stopper of the vial. Aseptically add one vial of the solution to every 200mL of modified sorbitol MacConkey (CT-SMAC) agar medium.

Item	Specification
Packing Size:	1 mL x 10 Vials/box
Formulation:	10 µg Cefixime, 0.5 mg Potassium Tellurite
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.