
Product Specification

Product Name: Nalidixic acid
Product Code: SG7333 / SG7279 / SG7276

PRODUCT DESCRIPTION:

Acts as a supplement for Tryptone Soy Broth and Listeria Enrichment Broth Base.

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial of 2.25 mg (SG7333) and one vial of 2.25 mg Polymyxin B (SG7332) to every 225mL of Tryptone Soy Broth medium to prepare Tryptone Soy Broth. Add one vial each of 4 mg (SG7279) and 5 mg Acriflavine (SG7278) to every 200 mL of Listeria Enrichment Broth Base medium to prepare LB2 medium. Add one vial each of 5 mg (SG7276) and 3 mg Acriflavine (SG7277) to every 225mL of Listeria Enrichment Broth Base medium to prepare LB1 medium.

Item	Specification
Packing Size:	1 mL x 10 Vials/box
Formulation:	5 mg / 4 mg / 2.25 mg Naphthylpyridinone Acid per mL
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be treated by autoclaving at 121°C for 30 minutes.