
Product Specification

Product Name: PALCAM Selective Additives

Product Code: SG7282

PRODUCT DESCRIPTION:

Supplement for PALCAM Agar Base, containing polymyxin B sulfate, acriflavine hydrochloride and ceftazidime.

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial of the solution to every 100 mL of PALCAM Agar Base medium.

Item	Specification
Packing Size:	1 mL x 10 Vials/box
Formulation:	1.0mg Polymyxin B sulfate, 0.5mg Acriflavine hydrochloride, 2.0mg Ceftazidime
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.