

---

## Product Specification

---

**Product Name:** Polymyxin B  
**Product Code:** SG7352 / SG7284 / SG7370 / SG7332

---

### PRODUCT DESCRIPTION:

Acts as a supplement for PALCAM, MYP, SPS Agar Bases, and Tryptone Soy Broth. Polymyxin B inhibits other accompanying gram-negative bacteria in the respective culture media.

### USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial of 1 mg (SG7352), 1 mg (SG7284), 1.2 mg (SG7370) or 2.25 mg (SG7332) to every 95 mL of MYP Agar Base, 100 mL of PALCAM Agar Base, 100 mL of SPS Agar Base or 225 mL of Tryptone Soy Broth respectively to prepare the respective medium.

---

Item	Specification
<b>Packing Size:</b>	1 mL x 10 Vials/box
<b>Formulation:</b>	1.0 mg (SG7352, SG7284)/1.2 mg (SG7370)/2.25 mg (SG7332) Polymyxin B
<b>Storage Condition:</b>	2-8°C
<b>Shelf Life:</b>	1 Year

---

**Waste Disposal:** After usage, bacteria-contaminated items should be autoclaved at 121°C for 30 minutes.