
Product Specification

Product Name: Rappaport Vassiliadis Soya Broth (RVS Broth)

Product Code: SG7261 / SG7262

PRODUCT DESCRIPTION:

Rappaport Vassiliadis Soya Broth (RVS Broth) is used for the selective enrichment of Salmonella spp. This product is available in both dehydrated and granulated form.

Use:

Stir to suspend 26.6 g of dehydrated/granulated medium in 1L distilled H₂O. Mix well with heat with frequent agitation. Boil for one minute. Autoclave for 15 minutes at 121°C to sterilize. Allow to cool (40 – 45°C) before making additions, such as antibiotics (if desired), before dispensing. For laboratory use only.

Components:

Soy Peptone – 4.5g/L

Sodium Chloride – 7.2g/L

Potassium Dihydrogen Phosphate – 1.26g/L

Dipotassium Phosphate – 0.18g/L

Anhydrous Magnesium Chloride – 13.4g/L

Malachite Green – 0.036g/L

Item

Description

Physical and Chemical Tests

Appearance: Yellowish Powder / Granules

pH (@25°C): 5.0 – 5.4

Storage: Store in a dry and dark place at room temperature. Tighten the cap immediately after use.