
Product Specification

Product Name: Salmonella Chromogenic Agar (SA)

Product Code: SG7140

PRODUCT DESCRIPTION:

Salmonella Chromogenic Agar (SA) is used for the detection and isolation of Salmonella species, including S. Typhi and S. Paratyphi.

Use:

Stir to suspend 34.9 g of dehydrated medium in 1L of distilled or deionized Water. Mix well with heat with frequent agitation. Boil for one minute until the powder is completely dissolved. Do not overheat or autoclave. For laboratory use only.

Components:

Peptone and yeast extract – 7.0g/L

Chromogenic and selective mix – 12.9g/L

Agar – 15.0g/L

Item	Description
------	-------------

Physical and Chemical Tests

Appearance: Yellowish Powder

pH (@25°C): 7.6±0.2

Storage: Store in a dry and dark place at room temperature. Tighten the cap immediately after use.