

---

## Product Specification

---

**Product Name:** Sulfadiazine Sodium

**Product Code:** SG7371

---

### PRODUCT DESCRIPTION:

Acts as a supplement for SPS Agar Base (Sulfite Polymyxin Sulfadiazine).

### USE:

Add one vial of the Sulfadiazine Sodium and one vial of 1.2 mg Polymyxin B (SG7370) to every 100 mL of SPS Agar Base medium.

---

---

Item	Specification
<b>Packing Size:</b>	1 mL x 10 Vials/box
<b>Formulation:</b>	12mg Sulfadiazine Sodium/vial
<b>Storage Condition:</b>	2-8°C
<b>Shelf Life:</b>	1 Year

---

**Waste Disposal:** After usage, bacteria-contaminated items should be treated by autoclaving at 121°C for 30 minutes.