
Product Specification

Product Name: Vitamin K1

Product Code: SG7339

PRODUCT DESCRIPTION:

Acts as a supplement for Peptone Yeast Glucose (PYG) Broth.

USE:

Before opening the ampoule, disinfect the surface with a 75% alcohol swab. Open the ampoule under sterile conditions. Add one vial each of the Vitamin K1 and Hemin (SG7338) to every 100 mL of PYG Base medium.

Item	Specification
Packing Size:	1 mL x 10 Vials/box
Formulation:	1.0 mg Vitamin K1 /mL
Storage Condition:	2-8°C
Shelf Life:	1 Year

Waste Disposal: After usage, bacteria-contaminated items should be treated by autoclaving at 121°C for 30 minutes.