
Product Specification

Product Name: Vibrio Chromogenic Agar

Product Code: SG7145

PRODUCT DESCRIPTION:

Vibrio Chromogenic Agar is used for isolating and identifying *Vibrio parahaemolyticus*, *Vibrio vulnificus*, and *Vibrio cholerae*.

Use:

Stir to suspend 74.7 g of dehydrated medium in 1L of distilled or deionized Water. Mix well with heat with frequent agitation. Boil to completely dissolve the powder, then cool to about 50°C. Do not overheat or autoclave. Pour into Petri dishes and set aside for use. For laboratory use only.

Components:

Peptone and yeast powder – 8.0g/L

Agar – 15.0g/L

Salts – 51.4g/L

Mixed colorants – 0.3g/L

Item	Description
------	-------------

Physical and Chemical Tests

Appearance: Yellowish Powder

pH (@25°C): 9.0±0.2

Storage: Store in a dry and dark place at room temperature. Tighten the cap immediately after use.