

Western Blot Filter Paper

Catalog #

CLW0003

Size

10 cm*15 cm

20 cm*15 cm

Description

Western Blot Filter Paper is used for the transfer of proteins from a gel onto a membrane

Components

100 pcs, Western Blot Filter Paper are made out of high-purity cotton fibers.

Thickness

180g/m²

Storage/Stability

Store at RT for 2 years.

COHESION BIOSCIENCES LIMITED

WEB

www.cohesionbio.com

ORDER

order@cohesionbio.com

SUPPORT

techsupport@cohesionbio.com

CUSTOM

custom@cohesionbio.com